

## Products for Biotechnology

### With Magnetic Porous Glass (MPG®)

**Protocol No.:** 5.105  
**Product:** MPG® Streptavidin (10 mg/ml, 4-6 × 10<sup>7</sup> particles/ml)  
**Procedure:** Purification of Nucleic Acids  
**Storage:** Refrigerate

PRODUCT NUMBER	DESCRIPTION	VOLUME
MSTR0502	MPG® Streptavidin, 5 µm, 50 nm (500 Å) pore diameter	2 ml (20 mg)
MSTR0510		10 ml (100 mg)

### General Procedure

**Materials:** (Based on 1.0 mg of MPG® Streptavidin, suspended in PBS, pH 7.5, 0.1% BSA, 0.02% NaN<sub>3</sub>)

DNA of Interest	EDTA, disodium (C <sub>10</sub> H <sub>14</sub> N <sub>2</sub> O <sub>8</sub> Na <sub>2</sub> )
Sodium Chloride (NaCl)	Low Speed Rotator
Potassium Chloride (KCl)	70°C Water Bath
Sodium Hydroxide (NaOH)	1.5 ml Nuclease-free Microcentrifuge Tubes
Sterile, Nuclease-Free, Deionized Water (dH <sub>2</sub> O)	Nuclease-free Pipettes and Pipette Tips
Tris Hydrochloride (Tris-HCl)	Magnetic Particle Separator, Prod. No. MPS0301 or MPS0001
Biotinylated Complementary DNA Probe	Vortex Mixer

#### Solution

2X Binding Buffer:  
(2 M KCl)

Washing Buffer  
(2 M NaCl)

2X Hybridization Buffer  
(1 M NaCl)

Storage Buffer  
(TE Buffer)

#### Preparation

Dissolve 15 g KCl in 80 ml dH<sub>2</sub>O. Adjust to pH 7.5 and bring volume to 100 ml with dH<sub>2</sub>O.

Dissolve 11.6 g NaCl in 80 ml dH<sub>2</sub>O. Adjust to pH 7.5 and bring volume to 100 ml with dH<sub>2</sub>O.

Dissolve 5.8 g NaCl in 80 ml dH<sub>2</sub>O. Adjust to pH 7.5 and bring volume to 100 ml with dH<sub>2</sub>O.

Dissolve 158 mg Tris-HCl and 37.2 mg EDTA in 95 ml of dH<sub>2</sub>O. Adjust to pH 7.5 with NaOH and bring volume to 100 ml with dH<sub>2</sub>O.

### Preparation of MPG® Streptavidin Complementary DNA Probe

- Transfer 100 µl (1.0 mg) of MPG® Streptavidin into a 1.5 ml nuclease-free microcentrifuge tube. Magnetically separate the MPG® Streptavidin from the solution by placing the tube in a Magnetic Particle Separator for at least 30 seconds. Carefully remove the supernatant with a pipette while the tube remains in the magnetic particle separator.
- Add 100 µl of 2X Binding Buffer and vortex to mix. Magnetically separate and carefully remove the supernatant.
- In a separate tube combine 500 pmol Biotinylated Complementary DNA Probe, 50 µl 2X Binding Buffer and dH<sub>2</sub>O to 100 µl total volume. Add this mixture to the MPG® Streptavidin. Vortex well and mix for 3 to 5 minutes at room temperature. Magnetically separate and carefully remove the supernatant.

4. Add 100  $\mu$ l of Washing Buffer to the Complementary DNA probe-bound MPG<sup>®</sup> Streptavidin particles and vortex to mix. Magnetically separate and carefully remove the supernatant.


### Hybridization of DNA

1. In a separate tube, combine 50  $\mu$ l 2X Hybridization Buffer, DNA of Interest and dH<sub>2</sub>O to 100  $\mu$ l total volume. Place the tube in a 70°C water bath for 2 to 3 minutes to denature the secondary structure. Add this mixture to the Complementary DNA Probe-bound MPG<sup>®</sup> Streptavidin particles and incubate 5 minutes at room temperature. Magnetically separate and carefully remove the supernatant.
2. Add 100  $\mu$ l of Washing Buffer and mix. Magnetically separate and carefully remove the supernatant.

### Recovery of Hybridized DNA

1. Add 100  $\mu$ l of dH<sub>2</sub>O or Storage Buffer to the dsDNA-bound MPG<sup>®</sup> Streptavidin particles. Place in a 70°C water bath for 2 to 3 minutes to release the purified ssDNA of Interest. Magnetically separate and quickly, but carefully, transfer the released ssDNA in the supernatant to a second tube. Repeat one more time. The ssDNA is now ready for further manipulation.
2. The ssDNA may be stored dry, frozen at -20°C in Storage Buffer or in ethanol at -20°C.

### Schematic Diagram Detailing Purification of Nucleic Acids



**FOR TECHNICAL SERVICE ON THIS OR ANY OTHER PureBiotech PRODUCT CALL 866-252-7771 or e-mail us at [info@purebiotechllc.com](mailto:info@purebiotechllc.com).**

*For in vitro research use only.*

MPG<sup>®</sup> is a registered trademark of Millipore Corporation. Magnetic Porous Glass and certain applications in which it is used are covered by U.S. Patent 5,601,979; 5,610,274 and 5,734,020 owned by Millipore Corporation and licensed by PureBiotech, LLC.