



## DNA CONTAMINANT REMOVAL SOLUTION

DNA-OFF™ is a non-alkaline, non-corrosive and non-carcinogenic ready-to-use cleaning solution that is highly active against DNA contamination on any laboratory surface.

### Directions for use:

DNA-OFF™ should be stored at room temperature; at colder temperatures a precipitate may form which is easily brought into solution at 37°C. Keep the cap tightly closed when not in use. Wear gloves, safety glasses, and take precautions to avoid inhaling aerosols when using this product.

**To decontaminate laboratory surfaces:** Apply DNA-OFF™ directly to the lab surface. After a few minutes wipe thoroughly with a paper towel, rinse with water and dry with a clean paper towel.

**To decontaminate laboratory apparatus:** Generously apply DNA-OFF™ to a paper towel and wipe all exposed surfaces of the apparatus thoroughly. After a few minutes rinse with water and dry with a clean paper towel. To clean small parts, briefly soak them in DNA-OFF™, rinse with water and dry.

**To decontaminate plastic and glass vessels:** Add ample DNA-OFF™ to enable coating the entire surface of the vessel by swirling or vortexing. After a few minutes discard the solution and rinse vessels thoroughly two times with distilled water.

**To decontaminate pipettors:** Following manufacturers instructions; remove the shaft from the pipettor and remove seals and gaskets from the shaft. Soak the shaft for a few minutes in DNA-OFF™, rinse the shaft thoroughly with water, let dry and reassemble.

### Caution:

May cause eye or skin irritation; may be harmful if inhaled or swallowed. Wear appropriate personal protective equipment to avoid exposure to this product. Consult the *Safety Data Sheet* for additional information. **FOR DNA DECONTAMINATION IN LABORATORY RESEARCH SETTING ONLY**

**GUARANTEED FOR ONE YEAR FROM PURCHASE.**