

Taq-*FORCE*[™] Hot DNA Polymerase

Product No.: HTAQ0250		HTAQ0500	
<i>System Includes:</i>		<i>System Includes:</i>	
1 x 250 U	Taq- <i>FORCE</i> [™] Hot DNA Polymerase (5 U/μl)	2 x 250 U	Taq- <i>FORCE</i> [™] Hot DNA Polymerase (5 U/μl)
1 x 1.2 ml	10X NH ₄ Reaction Buffer (pH 8.3)	2 x 1.2 ml	10X NH ₄ Reaction Buffer (pH 8.3)
1 x 1.2 ml	MgCl ₂	2 x 1.2 ml	MgCl ₂

Taq-*FORCE*[™] Hot is a heat-activated form of native Taq. The enzyme provides extremely high specificity and eliminates non-specific reaction products, to give higher yields than standard Taq DNA polymerase. The enzyme must be heat-activated to acquire catalytic activity, usually at 95°C for 7-10 minutes to give an optimal enzyme:template ratio for most reactions. Although usually unnecessary, the specificity can be maximized for a particular template by empirically determining the optimal pre-incubation time at 95°C.

Unit Definition: One unit is the amount of enzyme that will incorporate 10 nmoles of dNTPs into an acid-insoluble form in 30 minutes at 72°C under the following assay conditions: 25 mM TAPS, pH 9.3 (25°C); 50 mM KCl; 2 mM MgCl₂; 0.2 mM each dATP, dGTP, dTTP and 0.1 mM radiolabeled dCTP; 0.25 mg/ml activated salmon sperm DNA; 1 mM β-mercaptoethanol.

Storage Buffer: Enzyme is supplied in 20 mM Tris-HCl (pH 7.5), 100 mM NaCl, 0.1 mM EDTA, 2 mM DTT, 50% glycerol, 0.1% Tween-20.

Storage Conditions: -20°C. *DO NOT STORE IN A FROST-FREE FREEZER.*

Quality Control: Endonuclease, nickase, or exonuclease activities were not detectable after 8 hours incubation, respectively, of 1 μg of lambda, pBR322, or *Hind* III digested lambda DNA at 72°C in the presence of 5 units of Taq-*FORCE*[™] Hot DNA Polymerase.

10X NH₄ Reaction Buffer (pH 8.3): 160 mM (NH₄)₂SO₄, 670 mM Tris-HCl (pH 8.3 at 25°C), 0.1% Tween-20. The pH of this buffer is optimized for use with Taq-*FORCE*[™] Hot. Use of other buffers may reduce activity.

Mg⁺⁺ Stock Solution: 50 mM MgCl₂

Reaction Conditions: The optimal conditions (incubation time, temperatures, conc. of enzyme, template DNA, primers, MgCl₂) depend on the system and must be determined empirically. **IMPORTANT: Spin vials briefly before use.**

Component	Volume	Final Concentration
10X NH ₄ Reaction Buffer (pH 8.3)	5 μl	1X
DNTPs Pre-Mixed (<i>Cat. #DNTP10</i>)	4 μl	0.2 mM
MgCl ₂	variable	1.5 – 2.5 mM (typical)
Primer	variable	0.1 – 1.0 μM (each)
Taq- <i>FORCE</i> [™] Hot DNA Polymerase	variable	0.01 – 0.05 U/μl
Template DNA	variable	cosmid:0.1-5ng, genomic 0.1-0.5μg
Sterile H ₂ O	variable	-----
Final Volume	50 μl	-----

FOR RESEARCH USE ONLY

Note: Some applications in which this product can be used may be covered by patents issued and applicable in the United States and certain other countries. Purchase of this product does not convey a license to perform any patented process.